



व्हेटर्नरी डायग्नोस्टिक सेंटर
अचूक निष्कर्ष; प्राण्यांचा उत्कर्ष



LK Vet Care Pvt. Ltd.
॥ पशुसंपन्नतायर्म् तंत्रोपासना ॥



Only 15-20 min. A simple blood/serum test can be done directly in the cowshed

To check for pregnancy
28-Day Rapid Visual
Pregnancy Kit

Visual Pregnancy Kit Introduction:

Bovine Pregnancy Rapid Test Kit is a PAG (Pregnancy Assisted Glycoprotein) test, a blood sample based pregnancy test for cattle. The test can be used to detect the pregnancy in just 15-20 minutes at the farmers cow-shed by taking whole blood, serum or plasma samples from cow anytime after 28-days from breeding. It can help the farmers get proper guidance at their cow-shed itself without the need to transport blood sample to the nearest veterinary lab.

Features:



An accurate pregnancy test can be done anytime after the 28th day of breeding/inseminating the cow.



Shows results in just 15-20 minutes using whole blood, serum and plasma samples



Inspection can be done easily directly in cowshed without veterinary lab



A pregnancy test does not require any other equipment



A reliable kit designed to determine pregnancy in cattle



Easy-to-use, cost-effective, and gives fast results



The result is interpreted based on the color-coded lines of the test strip



Safer than rectal palpation and less invasive than ultrasound devices

Tap to Visit :

@lkvetcare.in

@info@lkvetcare.in

f pashulab

Advantages :



Improved Fertility:

Early detection of pregnancy reduces the chances of unwanted pregnancies and improves the efficiency of the reproductive system.



Non-Pregnant:

Identifies cows that are not pregnant and the duration of milk yield potential can be increased.



Saves Money:

Helps to save on feed, healthcare and labor costs associated with management of non-pregnant cattle.



Genetic improvement:

Early pregnancy detection can enable fertility decisions. Moreover, selection for genetic traits increases the long-term value of the herd.

Specifications and Usage:

Methodology	Pregnancy Assisted Glycoprotein (PAG Based)
shelf life	12 Months
sensitivity	99.00%
antigen specificity	98.07%

Applicable for	Bovine
Sample type	whole blood, serum
storage	4-28 Deg cent
packing	Single Use Test Kit (Minimum 5 Kits)

How to take a sample for testing:



1
Take 2-3 ml of blood from the tail



2
Add 3 drops of blood to the cassette



Add 2 drops of diluent to the sample



Results in just 15-20 minutes

Test Interpretation:



Positive



Negative



Invalid



New Cassette

Contents of the test kit:

Rapid Visual Pregnancy Kit 1) Test Cassette 2) Dropper 3) Blood Collection Tube 4) Blood Collection Needle 5) EDTA Vial 6) Diluent



For Purchase of Brucella Test Kit Contact: LK Vet Care Pvt. Ltd., 3A, Rachna Avenue, FC Road, Deccan Gymkhana, Pune - 9552524949 **Also available at: Pandharpur Pashulab:** Wakratund Building, Near KBP College Chowk, Wakhri Road, Pandharpur - 8956210651